

BioDAQ[®]

Rats and Mice

Automated Episodic Food & Liquid Intake Monitor



Home Cage

The BioDAQ Food/Liquid Intake Monitor automatically tracks natural feeding and drinking behavior of single-housed animals in a familiar home cage environment. Animals can live chronically in this environment with a variety of enrichment. Cages and hardware are easily cleaned with standard facility procedures.

Reduced Stress

Feeding and drinking modules attach outside the home cage for easy access by both animals and researchers. When networked, real-time data can be accessed and monitored remotely, reducing the need for frequent study room disturbances.

Electronic Sensors

Feeding and drinking modules are equipped with an electronic sensor that monitors hopper weight second by second. The software detects stable weight changes that identifies each interaction with the food or liquid source as a feeding or drinking bout. An integrated environmental monitor automatically records light cycles, temperature, and humidity.

Modular Hardware Design

Modular hardware design allows the researcher to configure the animal's cage to their specific study design.

- Food Intake
- Liquid Intake
- Food Choice and Place Preference
- Liquid Choice and Taste Preference
- Accurate intake tracking for food/drug mixtures
- Compatible with Telemetry devices
- Optogenetic/Tethered Cage options
- Body Weight
- Intracerebral cannula-compatible hoppers
- Easy system expansion
- Low vibration rack - accurate measurements


BioDAQ[®] INTAKE MONITOR

GEN3 Data Analysis Software

Gate Scheduler Tool

Manage individual access


- Supports multiple gate programs simultaneously
- Allows users to schedule future access dates
- Manage individual access for each food and liquid hopper



User-Friendly Software

New Rapid Interface


- Easy Data Analysis
- Enables seamless import of individual animal data
- Streamlines sorting and analysis for greater efficiency
- Several systems can be run together
- Multiple studies can be run on one rack
- Data can be securely accessed remotely by numerous users at a time



Graphing Tool

Visual data

- Features an intuitive interface
- Vivid graphing visualizations
- Seamless data export for effortless analysis




HOME CAGE

phenotype assessment

BioDAQ GEN3 Electronic

Automated Episodic Intake (same rack and hardware)



MOUSE CAGES

Modified standard mouse cages and stands are used for singly housed mouse studies.



RAT CAGES

Modified standard rat cages are used for singly housed rat studies.



NO LEAKS

The BioDAQ liquid intake bottle is designed not to leak. The drinking bout data is the animal's fluid intake, not a combination of intake and dripping.




REDUCE SPILLAGE

BioDAQ Extended mouse food hoppers retain the crumbs so that the weight is recorded as uneaten.



Biological Data Acquisition



Automated Gate Controller

Program access by time or amount consumed



The BioDAQ Automated Gate Controller lets researchers schedule openings based on time or a specified food/liquid intake. Mounted to the sensor, it integrates with the manual gate mechanism for seamless operation.

- Automated control enhances accuracy
- Reduces stress
- Minimizes researcher visits
- Pair feeding
- Intermittent fasting
- Binge eating
- Remotely fast animals in dark cycle

Body Weight Scale


Our new body weight scale enables automated body weight measurements at regular intervals throughout the day. The enclosed, dimly lit space naturally attracts mice, encouraging them to enter the hopper with ease.

- Allows for frequent body weight throughout the day to correlate with food and water intake data
- Easily check animal health remotely
- When the cage is open mice will enter the weight hopper reducing animal stress and allowing the researcher to remove for handling, dosing or visual inspection



Unplugged

Manual Periodic Intake




BioDAQ Unplugged includes cages and hardware with an adapter instead of an electronic sensor.

- Manually measure periodic intake
- Reduce spillage
- Reduce researcher/animal interaction
- Limit stress

Easily upgraded to BioDAQ GEN3 Electronic capabilities.


BioDAQ home cages can be customized with one or multiple feeding and drinking modules, adapting to any study design. The attachable hardware accommodates both food hoppers and liquid bottles, offering flexibility for researchers. Each cage includes a standard top and water bottle.



Food



Body Weight




Food & Liquid




Slotted Lid




Food and Liquid CHOICE



Split Cage and Lid



Automated Gate




Telemetry Cages



Optogenetic Cage

Applications

Biological Data Acquisition



Publications

researchdiets.com/biodaq

1 Byun S, Maric I, Stina Borchers, Sotzen MR, Olekanma D, Hayes MR, et al. **From the pancreas to the amygdala: new brain area critical for ingestive and motivated behavior control exerted by amylin.** *iScience*. 2025 Feb 1;112040–0.

Intracerebral Cannula

2 Sanchez-Navarro MJ, Borner T, Reiner BC, Crist RC, Samson WK, Yosten GLC, et al. **GPR-160 Receptor Signaling in the Dorsal Vagal Complex of Male Rats Modulates Meal Microstructure and CART-Mediated Hypophagia.** *Nutrients*. 2023 May 11;15(10):2268.

3 Yin TC, Mittal A, Buscaglia P, Li W, Sebag JA. **Activation of amygdala prokineticin receptor 2 neurons drives the anorexigenic activity of the neuropeptide PK2.** *Journal of Biological Chemistry*. 2022 Dec 17;299(1):102814–4.

1 Douglass AM, Kucukdereli H, Madara JC, Wang D, Wu C, Lowenstein ED, et al. **Acute and circadian feedforward regulation of agouti-related peptide hunger neurons.** *Cell Metabolism*. 2024 Dec 1;37(3).

Tethered and Optogenetics

2 Cimino I, Kim H, Tung YCL, Pedersen K, Rimmington D, Tadross JA, et al. **Activation of the hypothalamic-pituitary-adrenal axis by exogenous and endogenous GDF15.** *Proceedings of the National Academy of Sciences of the United States of America*. 2021 Jul 6;118(27):e2106868118.

3 Tracy AL, Schurdak JD, Chambers JB, Benoit SC. **Aversion learning can reduce meal size without taste avoidance in rats.** *Obesity*. 2016 Feb 1;24(3):606–14.

1 Hou T, Su W, Chacon AN, Lin AH, Guo Z, Gong MC. **Feeding- and Light-Cycle Synergistically Regulate Mouse Blood Pressure Daily Rhythm via Bmal1-Dependent and Independent Mechanisms.** *Journal of Biological Rhythms*. 2025 Jan 8;40(1):76–90.

Telemetry

2 Hou T, Chacon AN, Su W, Katsumata Y, Guo Z, Gong MC. **Role of sympathetic pathway in light-phase time-restricted feeding-induced blood pressure circadian rhythm alteration.** *Frontiers in Nutrition*. 2022 Sep 8;9.

3 Jackson C, Peterson R, Compton D, Ulman E, Tiesma S, Brockway R. **Continuous glucose and food intake monitoring in the male ZSD rat: comparison of normal and high fat diets (1051.12).** *The FASEB Journal*. 2014 Apr 1;28(S1).

1 Francois M, Canal Delgado I, Shargorodsky N, Leu CS, Zeltser L. **Assessing the effects of stress on feeding behaviors in laboratory mice.** *eLife*. 2022 Feb 15;11.

Automated Gates

2 Hopp K, Catenacci VA, Dwivedi N, Kline TL, Wang W, You Z, et al. **Weight loss and cystic disease progression in autosomal dominant polycystic kidney disease.** *iScience*. 2021 Dec 1;25(1):103697–7.

3 Reichenbach A, Stark R, Méquinion M, Lockie SH, Lemus MB, Mynatt RL, et al. **Carnitine acetyltransferase (Crat) in hunger-sensing AgRP neurons permits adaptation to calorie restriction.** *The FASEB Journal*. 2018 Jun 22;32(12):6923–33.

1 Elfers CT, Blevins JE, Salameh TS, Lawson EA, Silva D, Kiselyov A, et al. **Novel Long-Acting Oxytocin Analog with Increased Efficacy in Reducing Food Intake and Body Weight.** *International Journal of Molecular Sciences*. 2022 Sep 24;23(19):11249.

Meal Pattern Analysis

2 Yang TY, Gardner JC, Gentile JD, Liang NC. **Sex and individual differences in meal patterns mediate the persistency of running-associated high-fat diet avoidance in rats.** *American journal of physiology Regulatory, integrative and comparative physiology* 2019 Jan;316(2):R130–43.

3 Wheeler EC, Choi P, De Howitt J, Gill S, Watson S, Yu S, et al. **Cannabis Sativa targets mediobasal hypothalamic neurons to stimulate appetite.** *Scientific Reports [Internet]*. 2023 Dec 27;13(1):22970.

1 Pietrucci CL, Milton LK, Greaves E, Stefanidis A, van den Buuse M, Oldfield BJ, et al. **The BDNF Val66Met Polymorphism Does Not Increase Susceptibility to Activity-Based Anorexia in Rats.** *Biology*. 2022 May 1;11(5):623.

Preference Testing

2 Talukdar S, Owen Bryn M, Song P, Hernandez G, Zhang Y, Zhou Y, et al. **FGF21 Regulates Sweet and Alcohol Preference.** *Cell Metabolism*. 2016 Feb;23(2):344–9.

3 Chen YW, Barson JR, Chen A, Hoebel BG, Leibowitz SF. **Hypothalamic peptides controlling alcohol intake: Differential effects on microstructure of drinking bouts.** *Alcohol*. 2014 Nov;48(7):657–64.

RESEARCH DIETS INC.®

20 Jules Lane | New Brunswick, NJ 08901 USA | Tel: 732.247.2390 Fax: 732.247.2340
www.researchdiets.com/biodaq | info@researchdiets.com

© 2025 Research Diets, Inc. All rights reserved. BioDAQ GEN3-1500-10-25

